

Waterproof Digital Veterinary Thermometer

Hobday's Specification	°C	Man	Cat	Horse	Dog	Cow	Pig	Sheep	Fowl
		37	38	38	38.8	38.9	39.5	40	41.7

Instructions for use:

1. Press the ON/OFF button next to the display; a short alarm will sound, indicating that the thermometer is operational.
2. When the letters "Lo" and a flashing "C" display, the thermometer is ready for use.
3. Position the probe.
4. During the reading, the current temperature is displayed continuously and the "°C" symbol flashes.
5. When a constant temperature is reached, the measurement is completed. The highest temperature measured appears on the LCD, an alarm will sound four times and the "°C" symbol will stop flashing. The temperature is considered constant when the temperature rises less than 0.1°C within 4 seconds.
6. The unit will automatically turn off in 8-10 minutes. However, to prolong battery life, press the power switch once the temperature had been noted to turn off thermometer.
7. To switch between °C & °F, when turning the thermometer on, hold the ON/OFF button for 2 seconds.

Battery replacement:

The battery is empty and needs replacing when the "  " or "  " battery symbol appears on the right of the LCD. Remove the battery cover and replace it with a battery (preferably non-mercury) of the same type.

Please note:

The "+" sign up and "-" sign down.

Product disposal

Batteries do not belong in the domestic waste. Please hand them in at a collection point or recycling centre as special waste.

Warning: Do not attempt to recharge the battery. Do not dispose of battery in or near fire as it may explode.

Swallowing a button battery may be fatal.

Waterproof Digital Veterinary Thermometer

Features:

- Waterproof veterinary thermometer
- Easy to read, digital liquid crystal display
- Compact, accurate and durable large scale integration device
- Very sensitive device, quick to obtain results
- ABS resin body presents no hazard compared to broken mercury glass thermometer

Specifications:

Range:	32 to 42.9°C (89.6°F to 109.2°F)
Accuracy:	+/- 0.1°C between 35.5 and 39°C, otherwise +/-0.2°C
Display:	LCD 3.5 digits
Battery:	1pc x LR41 button battery (or equivalent)
Battery Life:	100 hours of continuous operation

Precaution

1. Do not drop - unit is not shock proof
2. Do not soak the display in water - only the probe is waterproof
3. Do not boil - clean the unit by wiping with a dry cloth and sterilize the probe with ethyl (rubbing) alcohol
4. Do not allow the unit to come into contact with thinners or chemical solvents.
5. Do not store in direct sunlight or at high temperature
6. Do not bend the unit
7. Do not attempt to disassemble the unit except to replace the battery.