



Product Safety Data Sheet

Brannan PSDS: 006.0
Issue Date: Feb. 1999
Revision: 0 Page: 1 of 1

Product Name: LoTox™ Laboratory Thermometer

The liquid in glass thermometers produced by S. Brannan & Sons Limited, contain expansion media comprising substances and preparations which present no risk to health, safety or environment when encased in the glass bulb and capillary.

If the instrument is broken during use, minimal risk is presented because of the small volume of materials being used. The following information is supplied as an aid to further reduce risk. As the Company may change constituent materials, users should check that they are in possession of the most recent revision of this Safety Data Sheet.

1. Manufacturers Information:

Address:

S. Brannan & Sons Limited
Leconfield Industrial Estate
Cleator Moor
Cumbria CA25 5QE
United Kingdom

Telephone:

UK: 01946 816600
International: 0044 1946 816600
Out of office hours: +44 778 979 1777

2. Medium Ingredients

This product contains ingredients classified as hazardous. However, they are not present in sufficient quantity to warrant classifying the product hazardous.

Expansion Medium - Maximum Volume 0.62 ml. per instrument
Colouring Dye - Average content 1.55 mg per instrument
Flammability - > 55° C

3. First Aid Data

Skin/eye contact - wash with water. Seek medical attention.

4. Spillage

Absorb in sand. Transfer to safe open area for atmospheric evaporation.

5. Other Materials

Organic fillings may react violently with strong acids, alkalis and oxidising agents.



Brannan Thermometers & Instrumentation

Leconfield Industrial Estate • Cleator Moor • Cumbria • CA25 5QE • UK

tel: +44 (0)1946 816600 • fax: +44 (0)1946 816625 • email: sales@brannan.co.uk • web: www.brannan.co.uk

