



# Product Safety Data Sheet

**Product Name: Maximum Reading Freezer Thermometer**

Brannan PSDS: 007.0  
Issue Date: March 2005  
Revision: 0 Page: 1 of 1

The liquid in glass thermometers produced by S. Brannan & Sons Limited, contain expansion media comprising substances and preparations which present no risk to health, safety or environment when encased in the glass bulb and capillary.

If the instrument is broken during use, minimal risk is presented because of the small volume of materials being used. The following information is supplied as an aid to further reduce risk. As the Company may change constituent materials, users should check that they are in possession of the most recent revision of this Safety Data Sheet.

## 1. Manufacturers Information:

**Address:**

S. Brannan & Sons Limited  
Leconfield Industrial Estate  
Cleator Moor  
Cumbria CA25 5QE  
United Kingdom

**Telephone:**

UK: 01946 816600  
International: 0044 1946 816600  
Out of office hours: +44 778 979 1777

## 2. Medium Constituents

Characteristic	FZ1 Solution	Methylated Spirit in Aqueous Solution	Petroleum Fraction 120
Flammability	N/A	F.P. (+12°C)	F.P. (+2°C)
Selected Safety Phrases	(S24) Avoid contact with skin.	(S24) Avoid contact with skin.	(S24) Avoid contact with skin. (S62) If swallowed, do not induce vomiting. Seek medical advice.
Selected Risk Phrases	(R22) Harmful if swallowed.	(R20) Harmful by inhalation. (R22) Harmful if swallowed.	(R65) May cause lung damage if swallowed. (R11) Highly flammable
Maximum Volume Per Instrument	0.4ml	0.8ml	0.05ml
Dye Characterisation	Eurocert Greens (EC No. E142)		

## 3. First Aid Data

If any of the media is ingested, seek medical attention.  
Skin/eye contact - wash with water and seek medical attention.

## 4. Spillage

Disposal, consistent with local authority regulations.



### Brannan Thermometers & Instrumentation

Leconfield Industrial Estate • Cleator Moor • Cumbria • CA25 5QE • UK

tel: +44 (0)1946 816600 • fax: +44 (0)1946 816625 • email: sales@brannan.co.uk • web: www.brannan.co.uk

