

Brannan Instruction Sheet



Countdown & Count Up Digital Timer (28/213/0)

Countdown Timer

1. Reset timer by pressing "M" and "S" at the same time.
2. Press "M" to set minute timer.
3. Press "S" to set second timer
4. Press "START/PAUSE" to start or pause timer.

Count Up Timer

1. Reset timer by pressing "M" and "S" at the same time.
2. Press "START/PAUSE" to start count up.
3. Press "START/PAUSE" to stop count up.

Battery Installation

Remove battery cover on back of timer and insert 1 x LR44 battery.

Battery Warning - Remove battery if timer will not be used for an extended period.

Keep battery away from children.

Care of your Product

1. Avoid exposing timer to extreme temperature, water or severe shock.
2. Avoid contact with any corrosive materials such as alcohol or cleaning agents.



Please dispose of electrical devices in accordance with local regulation

www.brannan.co.uk

Rev.1.2021

Measurement for Life.

