



Product Data Sheet

Body Temperature Infrared Thermometer

Non-contact and hygienic infrared thermometer specially designed to take the body temperature of a person regardless of room temperature. The thermometer has a fast response time allowing a greater number of measurements to be taken in a short time.

Designed to eliminate the chance of cross-contamination, it is ideally suited for use in schools, businesses or other environments where a large number of people are coming into contact with one another.

Features

- Easy and safe
- Selectable °C or °F
- Set alarm value
- Selectable body or surface temperature
- Memorization of the last 32 measurements
- Supplied with battery, full instructions and protective bag
- Laser targeting
- Trigger lock
- Automatic data hold
- Backlit LCD display



Product Specifications

Measuring Range:	In Body Mode: 32 to 42.5°C (90 to 108°F) In Surface Temp Mode: 0 to 60°C (32 to 140°F)
Operating Temp:	10 to 40°C (50 to 104°F)
Storage Temp:	0 to 50°C (32 to 122°F)
Accuracy:	+/-0.3°C (or 0.54°F)
Response Time:	0.5 second
Relative Humidity:	<85%
Emissivity:	0.95 fixed value
Measuring Distance:	5 - 15cm
Resolution:	0.1°C (0.1°F)
Battery:	2 x AA
Power Off:	Automatic shut off after 7 seconds, approx
CE & ROHS/WEEE:	Approved
Hazard Information (SDS):	See www.brannan.co.uk for information
LCD Dimensions:	28 x 24mm
Product Weight:	Packaged - 328g Unpackaged - 221g (incl. batteries, bag & instructions)
Product Dimensions:	Packaged - 186 x 32 x 70mm Unpackaged - 149 x 77 x 43mm
Display Box Weight:	88g
Display Box Dimensions:	200 x 111 x 60mm
Plastic Insert Tray Weight:	19g

Description	Barcode	Product No
Body Temperature Infrared Thermometer	5 011405 114006	11/400/3

Please note: all information on this data sheet is correct at the time of printing, but our policy is one of continuous product development. We reserve the right to change specification without notice.

PDS 0105:16:R2

Brannan Thermometers & Instrumentation

Leconfield Industrial Estate • Cleator Moor • Cumbria • CA25 5QE • UK

tel: +44 (0)1946 816600 fax: +44 (0)1946 816625 email: sales@brannan.co.uk web: www.brannan.co.uk



0582