

# Product data sheet

✉ [sales@brannan.co.uk](mailto:sales@brannan.co.uk)

☎ +44 (0)1946 816600

🌐 [brannan.co.uk](http://brannan.co.uk)

## ASTMH003 - ASTM SOIL HYDROMETERS

The American Society for Testing and Materials (ASTM) Specification E 100 describes hydrometers using the API (American Petroleum Institute) scale and others in degrees Specific Gravity.

All hydrometers can be supplied with UKAS certification, certificate traceable to Recognised National Standard or Statement of Conformity, if required.

### Product Features:

- Adjusted at 68°F
- 280mm long
- Made in UK
- Supplied in protective, resealable plastic tube

### Product Specifications:

<b>Material:</b>	Soda glass Steel shot ballast Synthetic wax resin
<b>Ranges:</b>	See table
<b>Overall Length:</b>	280mm
<b>Bulb Diameter:</b>	30mm
<b>Approx Weight:</b>	111g (incl. 45g protective plastic tube)

ASTM number	Range °SG	Division	Product number
151H	0.995 to 1.038°SG	0.001	09/151/0

  

ASTM number	Range g/ml	Division	Product number
152H	-5 to 60g/ml	1.0	09/152/0



View me on the website

