

# Product data sheet

✉ [sales@brannan.co.uk](mailto:sales@brannan.co.uk)

☎ +44 (0)1946 816600

🌐 [brannan.co.uk](http://brannan.co.uk)

## SPH016 - ALCOHOL HYDROMETER - GAY LUSSAC HYDROMETER (ALCOHOLOMETER)

250mm long, measuring alcohol % by volume and adjusted at 15°C (60°F).  
All hydrometers can be supplied with UKAS certification, certificate traceable to  
Recognised National Standard or Statement of Conformity, if required.

### Product features:

- Measures alcohol by volume
- 250mm long
- Adjusted at 15°C (60°F)
- Made in UK
- Supplied in protective, resealable plastic tube

### Product specifications:

<b>Material:</b>	Soda glass, steel shot ballast, synthetic wax resin
<b>Range:</b>	0 to 100
<b>Divisions:</b>	1.0
<b>Bulb diameter:</b>	16mm
<b>Overall length:</b>	250mm
<b>Weight:</b>	29g (including protective plastic tube and instruction sheet)

Range	Divisions	Product no
0 to 100 % Alcohol	1.0	10/670/0



View me on  
the website

