


Product data sheet

 sales@brannan.co.uk

 +44 (0)1946 816600

 brannan.co.uk

SPH016 - ALCOHOL HYDROMETER - GAY LUSSAC HYDROMETER (ALCOHOLOMETER)

250mm long, measuring alcohol % by volume and adjusted at 15°C (60°F).
All hydrometers can be supplied with UKAS certification, certificate traceable to
Recognised National Standard or Statement of Conformity, if required.

Product features:

- Measures alcohol by volume
- 250mm long
- Adjusted at 15°C (60°F)
- Made in UK
- Supplied in protective, resealable plastic tube

Product specifications:

Material:	Soda glass, steel shot ballast, synthetic wax resin
Range:	0 to 100
Divisions:	1.0
Bulb diameter:	16mm
Overall length:	250mm
Weight:	29g (including protective plastic tube and instruction sheet)

Range	Divisions	Product no
0 to 100 % Alcohol	1.0	10/670/0



View me on
the website

