

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

SPH017 - ALCOHOL HYDROMETER - ALCOHOL SIKES HYDROMETER

255mm hydrometers for measuring ethanol/water mixtures using an arbitrary Sikes scale at 51°F. All hydrometers can be supplied with UKAS certification, certificate traceable to Recognised National Standard or Statement of Conformity, if required.

Product features:

- Measures ethanol/water mixtures
- Uses arbitrary Sikes scales
- Adjusted at 51°F
- Made in UK
- Supplied in protective, resealable plastic tube

Product specifications:

Material:	Soda glass, steel shot ballast, synthetic wax resin
Ranges:	See table
Max length:	255mm
Bulb diameter:	31mm
Accuracy:	+/- 0.2
Divisions:	0.2
Approx weight:	105g (including 45g protective plastic tube)

Range	Product no
0 to 20 A	10/810/0
0 to 20	10/811/0
20 to 40	10/812/0
40 to 60	10/813/0
60 to 80	10/814/0
80 to 100	10/815/0





View me on
the website



Product data sheet


 sales@brannan.co.uk


 [+44 \(0\)1946 816600](tel:+44(0)1946816600)

 brannan.co.uk

Product data sheet

 sales@brannan.co.uk

 [+44 \(0\)1946 816600](tel:+44(0)1946816600)

 brannan.co.uk