


Product data sheet

 sales@brannan.co.uk

 +44 (0)1946 816600

 brannan.co.uk

BSH007 - BS 1377 SOIL HYDROMETERS

This soil hydrometer is to BS 1377 for the determination of particle size distribution in soil. It is 345mm long and adjusted at 20°C. All hydrometers can be supplied with UKAS certification, certificate traceable to Recognised National Standard or Statement of Conformity, if required.

Product Features:

- Adjusted at 20°C
- 345mm long
- Made in UK
- Supplied in protective, resealable plastic tube

Product Specifications:

Material:	Soda glass Steel shot ballast Synthetic wax resin
Range:	0.995 to 1.030 g/ml
Adjusted:	20°C
Overall Length:	345mm
Bulb Diameter:	26mm
Approx Weight:	103g (incl. 47g protective plastic tube)

BS number	Range °SG	Divisions	Product number
1377	0.995 to 1.030°SG	0.001	10/681/0



View on website

