


Product data sheet

 sales@brannan.co.uk

 +44 (0)1946 816600

 brannan.co.uk

BSH007 - ASTM SOIL HYDROMETERS

This soil hydrometer is to BS 1377 for the determination of particle size distribution in soil. It is 345mm long and adjusted at 20°F. All hydrometers can be supplied with UKAS certification, certificate traceable to Recognised National Standard or Statement of Conformity, if required.

Product features:

- Adjusted at 20°F
- 345mm long
- Made in UK
- Supplied in protective, resealable plastic tube

Product specifications:

Material: Soda glass
Steel shot ballast
Synthetic wax resin

Range: 0.995 to 1.030 g/ml

Adjusted: 20°C

Overall length: 345mm

Bulb Diameter: 26mm

Approx weight: 103g (incl. 47g protective plastic tube)

BS Number	Range °SG	Division	Product number
1377	0.995 to 1.030°SG	0.001	10/681/0



View me on the website

