

BS 1377 Soil hydrometers



This soil hydrometer is to BS 1377 for the determination of particle size distribution in soil. It is 345mm long and adjusted at 20°C. All hydrometers can be supplied with UKAS certification, certificate traceable to Recognised National Standard or Statement of Conformity, if required.

Product features:

- Adjusted at 20°C
- 345mm long
- Made in UK
- Supplied in protective, resealable plastic tube

Standard product specifications:

Material	Soda glass, steel shot ballast, synthetic wax resin
Range	0.995 to 1.030 g/ml
Adjusted:	20°C
Overall length	345mm
Bulb diameter	26mm
Approx weight	103g (incl. 47g protective plastic tube)

Standard product numbers:

BS number	Range °SG	Divisions	Product number
1377	0.995 to 1.030°SG	0.001	10/681/0

View me on
the website



REV0.2026