

# Brannan Product Data Sheet

PDS0161.20.REV3

## Digital Ovulation Thermometer and Chart


Easy to use battery operated electronic thermometer complete with instructions and recording charts for daily readings.

### Product Features:

- Compact, accurate and durable
- Non hazardous
- Complete with battery, instructions and recording chart

### Product Specifications:

<b>Temperature Range:</b>	32°C to 42°C
<b>Ingress Protection:</b>	IP22
<b>Accuracy:</b>	+/-0.1°C between 35.5°C and 42°C, otherwise +/-0.2°C
<b>Resolution:</b>	0.1°C
<b>Approx Response Time*:</b>	Approx 60 seconds
<b>Expected Lifespan:</b>	Approx 5 years
<b>Battery:</b>	1 x LR41 or equivalent (included)
<b>Approx Battery Life:</b>	100 hours continuous use (2+ years when idle)
<b>Storage Temperature:</b>	-25°C to +55°C / <95% RH
<b>Operating Temperature:</b>	+5°C to +35°C / <80% RH
<b>CE &amp; RoHS/WEEE:</b>	Compliant
<b>Hazard Information (SDS):</b>	See <a href="http://www.brannan.co.uk">www.brannan.co.uk</a> for information
<b>Product Dimensions:</b>	Packaged: 330mm x 180mm x 20mm Unpackaged: 210mm x 150mm x 18mm
<b>LCD Dimensions:</b>	7mm x 20mm
<b>Net Weight:</b>	58g (including batteries and instruction sheet)
<b>Display Card Weight:</b>	8g
<b>Polythene Bag Weight:</b>	10g
<b>Gross Weight:</b>	75g



Description	Barcode	Product No
Ovulation Kit	5 011405 110206	11/020/3

[www.brannan.co.uk](http://www.brannan.co.uk)



**Measurement for Life.**

