

Brannan Product Data Sheet

PDS058.20.REV6



Clinical Thermometer

Easy to use waterproof, electronic thermometer suitable for all ages, especially babies, young children and the elderly. Protective probe covers available for use with this thermometer.

Product Specifications:

Temperature Range:	32°C to 42.9°C
Ingress Protection:	IP22
Accuracy:	+/-0.1°C between 35.5°C and 42°C, otherwise +/-0.2°C
Resolution:	0.1°C
Approx Response Time*:	Approx 60 seconds
Expected Lifespan:	Approx 5 years
Battery:	1 x LR41 or equivalent (included)
Approx Battery Life:	100 hours continuous use (2+ years when idle)
Storage Temperature:	-25°C to +55°C / <95% RH
Operating Temperature:	+5°C to +35°C / <80% RH
CE & RoHS/WEEE:	Compliant
Hazard Information (SDS):	See www.brannan.co.uk for information
Product Dimensions:	Packaged: 167mm x 46mm x 16mm Unpackaged: 133mm x 24mm x 14mm
LCD Dimensions:	7mm x 20mm
Component Weights:	Net: 22g (including batteries and probe case) Display Box & Instructions: 8g
Gross Weight:	30g

Description	Barcode	Product No
Clinical Thermometer	5 011405 110640	11/064/2
Pack of 25pcs Disposable Probe Covers	5 011405 110664	11/066/0

www.brannan.co.uk



Measurement for Life.

