

Brannan Product Data Sheet

PDS0080.20.REV6

Infrared Ear Thermometer

A quick and easy infrared thermometer suitable for all ages.
Electronic thermometer with digital display and memory.
Complete with protective bag and replacement probe covers.
Especially suited with use with babies, young children and the elderly.
Pack of 20 replacement probe covers also available.

Features:

- Easy and safe
- Selectable °C or °F
- Supplied with battery, full instructions and protective bag

Product Specifications:

Measuring Temperature Range:	+35.5 to +42°C&F
Displayed Temperature Range:	+34 to +44°C&F
Operating Temperature Range:	+16 to +35°C (+60.8 to +95.0°F)
Operating Relative Humidity:	<80%RH
Resolution:	0.1°C&F
Accuracy:	+/-0.2°C between 35 and 42°C (+/-0.4°F between 95.9 and 107.6°F)
Battery:	2 x AAA
CE, ROHS/WEEE:	Approved
Hazard Information (SDS):	See www.brannan.co.uk for information
Product Weight:	Packaged - 175g Thermometer - 67g (including batteries) Probe Covers - 18g Display Box Weight: 75g Product Dimensions: LCD - 10mm x 20mm Thermometer - 140mm x 70mm x 40mm Packaged Product - 173mm x 116mm x 49mm



Description	Barcode	Product No
Infrared Ear Thermometer	5 011405 114204	11/420/3
Pack of 20pcs Disposable Probe Covers	5 011405 114211	11/421/0

www.brannan.co.uk



Measurement for Life.

