

Brannan Product Data Sheet

PDS0170.20.REV4

Forehead Strip Temperature Indicator

Quick and safe, forehead strip thermometer using heat sensitive liquid encased in plastic film. Ideal for small children or for adding into first aid kit. Supplied complete with protective case and full instructions and available packed in individual packaging or bulk packed.

Product Features:


- Simple to use and easy to read
- Quick and safe
- Supplied with storage box
- Made in UK

Product Specifications:

Temperature Range:	+35 to +40°C&F
Indicators:	Every 1.0°C/°F
Accuracy:	+/- 1.0°C
Product Material:	Thermometer strip: Terylene (PET)
Storage & Operating Conditions:	16 to 40°C (60 to 104°F) at 15 to 95% RH*
CE & RoHS/WEEE:	Compliant
Hazard Information (SDS):	See www.brannan.co.uk for information
Product Dimensions:	Net: 39 x 99 x 10mm (case closed) Gross: 157 x 85 x 14mm (display packed)
Component Weights:	Forehead Strip: 0.5g Plastic Case: 10g Paper Instructions: 2g Paper Instructions & packaging: 9.5g Gross: 22g (display packed)

*Note:

After use, store thermometer in the protective case and keep away from Ultraviolet/sunlight, dampness, or heat. The thermometer can be cleaned by gently wiping with a soft cloth or tissue dipped in water. Do not immerse thermometer in liquids
Optimal storage and usage between 16 to 27°C (60 to 80°F) at 40 to 60% RH.



Description	Barcode	Product No
Forehead Temperature Indicator - Carded	5 011405 110305	11/465/3
Forehead Temperature Indicator - Bulk	-	11/466/3

www.brannan.co.uk



Measurement for Life.

