



# Classroom Kit

This kit has been developed for use in the classroom. It comes complete with a variety of laboratory thermometers and a set of 'Classroom' notes and photocopiable worksheets to assist teachers to operate within the National Curriculum.

**CE/ROHS/WEEE:** N/A  
**Hazard Information (SDS):** See [www.brannan.co.uk](http://www.brannan.co.uk)  
**Component Weights:** Glass - 296g  
 Paper - 136.5g  
 Plastic - 133.5g  
**Gross Weight:** 566.0g  
**Product Dimensions :** 300 x 215 x 75mm



## Contents

- 8 x 150mm *initial*<sup>TM</sup> brand green filled thermometers
  - 7 x 300mm *initial*<sup>TM</sup> brand green filled thermometers
  - 1 x 300mm uncalibrated green filled thermometers
- See following page for details of *initial*<sup>TM</sup> brand laboratory thermometers.

Description	Barcode	Product No
Classroom Kit	5 011405 257000	25/700/0



[www.brannan.co.uk](http://www.brannan.co.uk)

Leconfield Industrial Estate Cleator Moor Cumbria CA25 5QE UK

**Est. 1913**  
Over 100 years  
of Thermometer  
Expertise




# Classroom Kit



The **initial**<sup>TM</sup> Brand range of laboratory thermometers offer a significant saving on our standard range of thermometers while making no compromises to the traditional **Brannan** quality and accuracy. This entry level range is now available in a choice of green or blue spirit filling and have Tri-top, moulded glass non-roll device to help prevent breakage. These thermometers are printed with a clear, easy-to-read scale making them particularly suitable for student use.



The **Brannan breaksafe**<sup>TM</sup> range of laboratory thermometers are encapsulated in PVC to contain any glass or filling material in the event of breakage.



Range °C only	Division	Length	Filling Colour	76mm Immersion		Total Immersion	
				Accuracy	Product No	Accuracy	Product No
-10 to +110°C	1.0°C	150mm	Green	-	-	+/-2.0°C @ 35°C	46/330/5
-10 to +110°C	1.0°C	150mm	Blue	-	-	+/-2.0°C @ 35°C	46/330/7
-10 to +110°C	1.0°C	200mm	Green	+/-1.5°C @ 35°C	46/323/5	-	-
-10 to +110°C	1.0°C	200mm	Blue	+/-1.5°C @ 35°C	46/323/7	-	-
-10 to +110°C	1.0°C	300mm	Green	+/-1.5°C @ 35°C	46/303/5	-	-
-10 to +110°C	1.0°C	300mm	Blue	+/-1.5°C @ 35°C	46/303/7	-	-
<b>PVC Coated <i>initial</i><sup>TM</sup> Brand* Thermometers - Green Filling*</b>							
-10 to +110°C	1.0°C	200mm	Green	+/-1.5°C @ 35°C	44/425/0	-	-
<b>Uncalibrated Thermometers (for student calibration)</b>							
Covers Range -10 to +110°C		305mm	Green	-	46/603/5	-	-
		305mm	Blue	-	46/603/7	-	-

\* For details about the full range of PVC coated thermometers, please see previous page of this catalogue