

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

CLASSROOM KIT - LABORATORY THERMOMETER KIT

This kit has been developed for use in the classroom. It comes complete with a variety of laboratory thermometers and a set of 'Classroom' notes and photocopiable worksheets to assist teachers to operate within the National Curriculum.

Product features:

- Selection of laboratory thermometer for pupils' investigative work
- Comprehensive worksheets for various experiments
- Made in UK

Product specifications:

CE/ROHS/WEEE:	N/A
Hazard information (SDS):	See brannan.co.uk
Component weights:	Glass: 296g Paper: 136.5g Plastic: 133.5g
Gross weight:	566.0g
Product dimensions :	300mm x 215mm x 75mm

Contents:

- 8 x 150mm *initial*TM brand green filled thermometers
 - 7 x 300mm *initial*TM brand green filled thermometers
 - 1 x 300mm uncalibrated green filled thermometers
- See following page for details of *initial*TM brand thermometers.



Description	Barcode	Product no
Classroom kit - laboratory thermometer kit	5 011405 257000	25/700/0

Product data sheet

✉ sales@brannan.co.uk


☎ +44 (0)1946 816600

🌐 brannan.co.uk

The *initial*TM Brand range of laboratory thermometers offer a significant saving on our standard range of thermometers while making no compromises to the traditional **Brannan** quality and accuracy. This entry level range is now available in a choice of green or blue spirit filling and have Tri-top, moulded glass non-roll device to help prevent breakage. These thermometers are printed with a clear, easy-to-read scale making them particularly suitable for student use.



The **Brannan breaksafe**TM range of laboratory thermometers are encapsulated in PVC to contain any glass or filling material in the event of breakage.



Range °C only	Division	Length	Filling Colour	76mm Immersion		Total Immersion	
				Accuracy	Product No	Accuracy	Product No
-10 to +110°C	1.0°C	150mm	Green	-	-	+/-2.0°C @ 35°C	46/330/5
-10 to +110°C	1.0°C	150mm	Blue	-	-	+/-2.0°C @ 35°C	46/330/7
-10 to +110°C	1.0°C	200mm	Green	+/-1.5°C @ 35°C	46/323/5	-	-
-10 to +110°C	1.0°C	200mm	Blue	+/-1.5°C @ 35°C	46/323/7	-	-
-10 to +110°C	1.0°C	300mm	Green	+/-1.5°C @ 35°C	46/303/5	-	-
-10 to +110°C	1.0°C	300mm	Blue	+/-1.5°C @ 35°C	46/303/7	-	-
PVC Coated <i>initial</i>TM Brand* Thermometers - Green Filling*							
-10 to +110°C	1.0°C	200mm	Green	+/-1.5°C @ 35°C	44/425/0	-	-
Uncalibrated Thermometers (for student calibration)							
Covers Range -10 to +110°C	305mm	Green	-	46/603/5	-	-	
	305mm	Blue	-	46/603/7	-	-	

* For details about the full range of PVC coated thermometers, please see previous page of this catalogue

View me on the website

