

# Product data sheet

✉ [sales@brannan.co.uk](mailto:sales@brannan.co.uk)

☎ +44 (0)1946 816600

🌐 [brannan.co.uk](http://brannan.co.uk)

## CLASSROOM KIT - LABORATORY THERMOMETER KIT



This kit has been developed for use in the classroom. It comes complete with a variety of laboratory thermometers and a set of 'Classroom' notes and photocopyable worksheets to assist teachers to operate within the National Curriculum.

### Product features:

- Selection of laboratory thermometer for pupils' investigative work
- Comprehensive worksheets for various experiments
- Made in UK

### Product specifications:

<b>CE/ROHS/WEEE:</b>	N/A
<b>Hazard information (SDS):</b>	See <a href="http://brannan.co.uk">brannan.co.uk</a>
<b>Component weights:</b>	Glass: 296g Paper: 136.5g Plastic: 133.5g
<b>Gross weight:</b>	566.0g
<b>Product dimensions :</b>	300mm x 215mm x 75mm

### Contents:

- 8 x 150mm *initial*<sup>™</sup> brand green filled thermometers
  - 7 x 300mm *initial*<sup>™</sup> brand green filled thermometers
  - 1 x 300mm uncalibrated green filled thermometers
- See following page for details of *initial*<sup>™</sup> brand thermometers.

Description	Barcode	Product no
Classroom kit - laboratory thermometer kit	5 011405 257000	25/700/0

# Product data sheet

✉ [sales@brannan.co.uk](mailto:sales@brannan.co.uk)


☎ +44 (0)1946 816600

🌐 [brannan.co.uk](http://brannan.co.uk)

The *initial*<sup>TM</sup> Brand range of laboratory thermometers offer a significant saving on our standard range of thermometers while making no compromises to the traditional **Brannan** quality and accuracy. This entry level range is now available in a choice of green or blue spirit filling and have Tri-top, moulded glass non-roll device to help prevent breakage. These thermometers are printed with a clear, easy-to-read scale making them particularly suitable for student use.



The **Brannan breaksafe**<sup>TM</sup> range of laboratory thermometers are encapsulated in PVC to contain any glass or filling material in the event of breakage.



Range °C only	Division	Length	Filling Colour	76mm Immersion		Total Immersion	
				Accuracy	Product No	Accuracy	Product No
-10 to +110°C	1.0°C	150mm	Green	-	-	+/-2.0°C @ 35°C	46/330/5
-10 to +110°C	1.0°C	150mm	Blue	-	-	+/-2.0°C @ 35°C	46/330/7
-10 to +110°C	1.0°C	200mm	Green	+/-1.5°C @ 35°C	46/323/5	-	-
-10 to +110°C	1.0°C	200mm	Blue	+/-1.5°C @ 35°C	46/323/7	-	-
-10 to +110°C	1.0°C	300mm	Green	+/-1.5°C @ 35°C	46/303/5	-	-
-10 to +110°C	1.0°C	300mm	Blue	+/-1.5°C @ 35°C	46/303/7	-	-
<b>PVC Coated <i>initial</i><sup>TM</sup> Brand* Thermometers - Green Filling*</b>							
-10 to +110°C	1.0°C	200mm	Green	+/-1.5°C @ 35°C	44/425/0	-	-
<b>Uncalibrated Thermometers (for student calibration)</b>							
Covers Range -10 to +110°C	305mm	Green	-	46/603/5	-	-	
	305mm	Blue	-	46/603/7	-	-	

\* For details about the full range of PVC coated thermometers, please see previous page of this catalogue

View me on the website

