

# Product data sheet

✉ [sales@brannan.co.uk](mailto:sales@brannan.co.uk)

☎ +44 (0)1946 816600

🌐 [brannan.co.uk](http://brannan.co.uk)

## FOREHEAD STRIP TEMPERATURE INDICATOR

Quick and safe, forehead strip thermometer using heat sensitive liquid encased in plastic film. Ideal for small children or for adding into first aid kit. Supplied complete with protective case and full instructions and available packed in individual packaging or bulk packed.

### Product features:


- Simple to use and easy to read
- Quick and safe
- Supplied with storage box
- Made in UK

### Product specifications:

|  |  |
|--|--|
| <b>Temperature range:</b>                  | +35 to +40°C&F   |
| <b>Indicators:</b>                         | Every 1.0°C/°F   |
| <b>Accuracy:</b>                           | +/- 1.0°C  |
| <b>Product material:</b>                   | Thermometer strip: Terylene (PET)  |
| <b>Storage &amp; operating conditions:</b> | 16 to 40°C (60 to 104°F) at 15 to 95% RH*  |
| <b>CE &amp; RoHS/WEEE:</b>                 | Compliant  |
| <b>Hazard Information (SDS):</b>           | See <a href="http://brannan.co.uk">brannan.co.uk</a> for information                             |
| <b>Product dimensions:</b>                 | Net: 39mm x 99mm x 10mm (case closed)<br>Gross: 157mm x 85mm x 14mm (display packed)             |
| <b>Product weights:</b>                    | Forehead strip: 0.5g<br>Plastic case: 10g<br>Paper instructions & packaging: 11.5g<br>Gross: 22g |

### \*Note:

After use, store thermometer in the protective case and keep away from Ultraviolet/sunlight, dampness, or heat. The thermometer can be cleaned by gently wiping with a soft cloth or tissue dipped in water. Do not immerse thermometer in liquids  
Optimal storage and usage between 16 to 27°C (60 to 80°F) at 40 to 60% RH.



| Description                             | Barcode         | Product No |
|---|-----------------|------------|
| Forehead Temperature Indicator - Carded | 5 011405 110305 | 11/465/3   |
| Forehead Temperature Indicator - Bulk   | -               | 11/466/3   |

[brannan.co.uk](http://brannan.co.uk)



View me on the website

