

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

ORAL CLINICAL THERMOMETER

Easy to use waterproof, electronic thermometer suitable for all ages, especially babies, young children and the elderly.
Protective probe covers available for use with this thermometer.

Product features:

- Compact, accurate and durable
- Large scale integration device
- Non hazardous
- Complete with battery and instructions
- Available pre-calibrated as part of **Brannan THERMcert™** range*



Product specifications:

Temperature range:	32°C to 42.9°C
Accuracy:	+/-0.1°C between 35.5°C and 42°C, otherwise +/-0.2°C
Resolution:	0.1°C
Approx response time*:	Approx 60 seconds
Expected lifespan:	Approx 5 years
Battery:	1 x LR41 or equivalent (included)
Approx battery life:	100 hours continuous use (2+ years when idle)
Storage temperature:	-25°C to +55°C / <95% RH
Operating temperature:	+5°C to +35°C / <80% RH
CE & RoHS/WEEE:	Compliant
LCD dimensions:	8mm x 18mm
Product dimensions:	Net: 133mm x 24mm x 14mm Gross: 167mm x 46mm x 16mm
Product weight:	Net: 22g (including batteries and instruction sheet)
Gross weight:	Display box and instructions: 8g Gross: 30g

Description	Barcode	Product No
Clinical Thermometer	5 011405 110640	11/064/2
Pack of 25pcs Disposable Probe Covers	5 011405 110664	11/066/0

Description	Calibration points	Product no
THERMcert™ pre-calibrated thermometer	Fixed at 37°C	TC-009*
Recalibration product reference: TC-009-RECALIBRATION		

* See website for details for Brannan **THERMcert™** pre-calibrated thermometers.



View me on
the website

