

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

ORAL CLINICAL THERMOMETER

Easy to use electronic thermometer suitable for all ages, especially babies, young children and the elderly. Protective probe covers available for use with this thermometer.

Product features:

- Compact, accurate and durable
- Large scale integration device
- Non hazardous
- Complete with battery and instructions
- Available pre-calibrated as part of **Brannan THERMcert™** range*



Product specifications:

Temperature range:	32°C to 42.9°C
Accuracy:	+/-0.1°C between 35.5°C and 42°C, otherwise +/-0.2°C
Resolution:	0.1°C
Approx response time*:	Approx 60 seconds
Expected lifespan:	Approx 5 years
Battery:	1 x LR41 or equivalent (included)
Approx battery life:	100 hours continuous use (2+ years when idle)
Storage temperature:	-25°C to +55°C / <95% RH
Operating temperature:	+5°C to +35°C / <80% RH
CE & RoHS/WEEE:	Compliant
LCD dimensions:	8mm x 18mm
Product dimensions:	Net: 133mm x 24mm x 14mm Gross: 167mm x 46mm x 16mm
Product weight:	Net: 22g (including batteries and instruction sheet)
Gross weight:	Display box and instructions: 8g Gross: 30g

Description	Barcode	Product No
Clinical Thermometer	5 011405 110640	11/064/2
Pack of 25pcs Disposable Probe Covers	5 011405 110664	11/066/0

Description	Calibration points	Product no
THERMcert™ pre-calibrated thermometer	Fixed at 37°C	TC-009*
Recalibration product reference: TC-009-RECALIBRATION		

* See website for details for Brannan **THERMcert™** pre-calibrated thermometers.



View me on the website

