

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

MINI DIAL VIVARIUM THERMOMETER


A compact instrument with a self adhesive pad which can be fixed inside a climate controlled chamber such as a vivarium.

Product features:

- 65mm diameter temperature meter
- Clear temperature scale on white moulding with red pointer
- Supplied with double sided adhesive pads for convenient fixture
- Made in UK

Product specifications:

Temperature range:	-10 to +50°C only or C&F
Divisions:	1°C / 2°F
Accuracy:	+/-2.0°C
CE & RoHS/WEEE:	N/A
Hazard information (SDS):	See brannan.co.uk for information
Product material:	Dial back: high impact styrene (6g) Dial window: clear crystal styrene (5.5g) Hairspring: bimetal sensor with copper knob (0.5g)
Product weight:	Net: 12g Display packaging: 24g Gross: 36g
Dimensions:	Net: ø65mm x 12mm Gross: 143mm x 110mm x 12mm



Description	Range	Barcode	Product No
Mini dial vivarium thermometer	-10 to 50°C only	5 011405 304407	30/440/2
	-10 to 50°C&F	5 011405 304407	30/440/3

View me on the website



brannan.co.uk

Measurement for life.

