

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

CYH001 - HYDROMETERS CYLINDERS

Soda glass hydrometer cylinders available with a standard internal diameter of 37mm and plastic foot.

Product features:


- Choice of depths for use with Brannan hydrometers
- Plastic foot for stability
- Made in UK
- Supplied in protective, resealable plastic tube

Product specifications:

Material: Soda glass
Plastic foot

See table for details

Internal Depth	Suitable for Hydrometer Series	Weight	Product Number
200mm	S50	88.5g	10/000/0
300mm	M50 & M100	144.0g	10/000/1
330mm	L50	157.5g	10/000/2



View me on the website

