Product data sheet





+44 (0)1946 816600



MASON'S WET & DRY BULB HYGROMETER

Product features:

 Simple psychrometer to measure water vapour in the atmosphere, on a percentage scale known as Relative Humidity (%RH)

 Includes table to calculate the difference between the "wet" and "dry" bulb temperature to give the %RH reading

Product specifications:

Temperature ranges: -5 to +50°C&F
Accuracy: +/- 1°C&F
Divisions: 1°C&F
CE/ROHS/WEEE: N/A

Hazard information (SDS): See brannan.co.uk for information

Product material: Case: polypropylene

Scale: polystyrene
Reservoir: polyethylene
Glass Tube: kerosene filled
Net: 310mm x 77mm x 32mm

Product dimensions: Net: 310mm x 77mm x 32mm

Gross: 315mm x 80mm x 32mm*

Product weight: Net: 110g

Box: 33g

Gross: 146g* (incl. chart)

^{*} Supplied in individual box

Description	Barcode	Product no
Masons hygrometer	5 011405 132222	13/222/3
Replacement reservoir pot	-	13/087/0



View me on the website



