

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

ADVANCED, NON-CONTACT INFRARED BODY & SURFACE THERMOMETER

Non-contact and hygienic infrared thermometer specially designed to take the body temperature of a person regardless of room temperature. The thermometer has a fast response time allowing a greater number of measurements to be taken in a short time.

Designed to eliminate the chance of cross-contamination, it is ideally suited for use in schools, businesses or other environments where a large number of people and coming into contact with one another.

Product features:

- Easy and safe
- Selectable °C or °F
- Selectable body or surface temperature
- Memorisation of the last 32 measurements
- Supplied with battery, full instructions and protective bag
- Trigger lock
- Automatic data hold
- LCD display
- Set alarm value

Product specifications:

Effective distance:	5 - 15cm
Body temperature range:	32°C to 42.5°C / 90 to 108°F
Surface temperature range:	0°C to 60°C / 32 to 140°F
Storage temperature:	0°C to +50°C&F
Operating temperature:	+10°C to +40°C&F
Accuracy:	+/-0.3°C / 0.54°F
Measuring time:	0.5 seconds
Relative humidity:	<85%
Emissivity:	0.95 fixed value
Resolution:	0.1°C
Battery:	2 x AA or equivalent (included) Voltage: 3.2V Work current: 8mA Continuous use time: approx 121 hours After approx 7 seconds inactivity
Auto power off:	Compliant
CE & RoHS/WEEE:	Compliant
LCD dimensions:	28mm x 24mm
Product dimensions:	Gross: 186mm x 132mm x 70mm Net: 149mm x 77mm x 43mm
Component weights:	Net: 221g (including batteries, bag & instructions) Display box: 88g Plastic insert tray: 19g Gross: 328g



Description	Barcode	Product No
Advanced body & surface, non-contact infrared thermometer	5 011405 114006	11/400/3

View me on
the website

