

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

DIGITAL OVULATION THERMOMETER AND CHART

Easy to use battery operated electronic thermometer complete with instructions and recording charts for daily readings.

Product features:

- Compact, accurate and durable
- Non hazardous
- Complete with battery, instructions and recording chart

Product specifications:

| | |
|---------------------------|--|
| Temperature range: | 32°C to 42°C |
| Ingress protection: | IP22 |
| Resolution: | 0.1°C |
| Approx response time*: | Approx 60 seconds |
| Accuracy: | +/-0.1°C between 35.5°C and 42°C, otherwise +/-0.2°C |
| LCD dimensions: | 7mm x 20mm |
| Expected lifespan: | Approx 5 years |
| Battery: | 1 x LR41 or equivalent (included) |
| Approx battery life: | 100 hours continuous use (2+ years when idle) |
| Storage temperature: | -25°C to +55°C / <95% RH |
| Operating temperature: | +5°C to +35°C / <80% RH |
| CE & RoHS/WEEE: | Compliant |
| Hazard information (SDS): | See brannan.co.uk for information |
| Product dimensions: | Net: 210mm x 150mm x 18mm Gross: 330mm x 180mm x 20mm |
| Product weights: | Net: 58g (including batteries and instruction sheet) Card: 8g Polythene bag: 10g Gross: 75g |



| Description | Barcode | Product No |
|---------------|-----------------|------------|
| Ovulation Kit | 5 011405 110206 | 11/020/3 |

View me on the website

