

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

DIGITAL OVULATION THERMOMETER AND CHART

Easy to use battery operated electronic thermometer complete with instructions and recording charts for daily readings.

Product features:

- Compact, accurate and durable
- Non hazardous
- Complete with battery, instructions and recording chart

Product specifications:

Temperature range:	32°C to 42°C
Ingress protection:	IP22
Resolution:	0.1°C
Approx response time*:	Approx 60 seconds
Accuracy:	+/-0.1°C between 35.5°C and 42°C, otherwise +/-0.2°C
LCD dimensions:	7mm x 20mm
Expected lifespan:	Approx 5 years
Battery:	1 x LR41 or equivalent (included)
Approx battery life:	100 hours continuous use (2+ years when idle)
Storage temperature:	-25°C to +55°C / <95% RH
Operating temperature:	+5°C to +35°C / <80% RH
CE & RoHS/WEEE:	Compliant
Hazard information (SDS):	See brannan.co.uk for information
Product dimensions:	Net: 210mm x 150mm x 18mm Gross: 330mm x 180mm x 20mm
Product weights:	Net: 58g (including batteries and instruction sheet) Card: 8g Polythene bag: 10g Gross: 75g



Description	Barcode	Product No
Ovulation Kit	5 011405 110206	11/020/3

View me on the website



brannan.co.uk

Measurement for life.

