

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

SPH004 - SPECIAL PURPOSE SPECIFIC GRAVITY HYDROMETER - SPECIFIC GRAVITY & QUEVENNE LACTOMETERS

These Lactometers are 160mm long milk hydrometers, adjusted at 15.6°C (60°F). Scales calibrated in specific gravity or Quevenne. 1° Quevenne equals 0.001°SG. All hydrometers can be supplied with UKAS certification, certificate traceable to Recognised National Standard or Statement of Conformity, if required.

Product features:

- Adjusted at 15.6°C (60°F)
- 160mm long
- Made in UK
- Supplied in protective, resealable plastic tube

Product specifications:

Material: Soda glass, steel shot ballast, synthetic wax resin
Ranges: See table below
Accuracy: See table below
Divisions: See table below
Overall length: 160mm
Bulb diameter: 16mm
Approx weight: 22g (including 8g protective plastic tube)

Range	Accuracy at 20°C	Division	Product number
1.025 to 1.035°SG	+/- 0.002	0.001	10/601/0
15 to 40° Quevenne/ 1.015 to 1.025°SG	+/- 2.000	1.000/ 0.100	10/660/0



View me on
the website

