

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

BENCH TOP TIMER WITH EXTERNAL TRIGGERING CONTACTS


The Brannan bench top timers are easy to use and are particularly suited for use in schools. They have colour coded buttons and a clear LCD display, designed to offer accurate measurement over a long period of time. This model has external triggering however timer is available without.

Product features:

- Three 4mm sockets for external triggering
- Large LCD and buttons
- Up to 9 hours, 59 minutes, 59 seconds
- 1/100 second resolution
- Includes batteries and operating instructions

Product specifications:

Timing period:	9 hours / 59 minutes / 59 seconds
Accuracy:	+/- 1 second
Resolution:	1/100 second
Product material:	ABS
Battery:	1pc x AA (supplied)
Operating temperature:	0 to 50°C
Storage temperature:	-10 to 60°C
CE & ROHS/WEEE:	Approved
Hazard information (SDS):	See brannan.co.uk for information
LCD dimensions:	63.7mm x 19.7mm
External contacts:	4mm sockets and 1m cables
Product dimensions:	Net: 100mm x 97mm x 70.5mm Gross: 138mm x 123mm x 80mm
Product weight:	Net: 257g Display box: 138g Gross: 395g



Description	Barcode	Product no.
Bench top timer with external triggering contacts	5 011405 284105	28/410/0
Spare contacts	n/a	28/410/4
Bench top timer without contacts	5 011405 284112	28/411/0



View me on the website



brannan.co.uk

Measurement for life.

