

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

SLIMLINE TIMER

Easy to use 99 minute / 59 second countdown and up timer with magnetic fixing.


Supplied with battery and full instructions in individual box.

Product features:

- Easy to use
- Count Up and down timers
- Clear audible alarm
- Magnetic fixture
- Battery and instructions included

Product specifications:

Timer count:	99 minute / 59 second
Product material:	ABS
Duration of alarm:	30 seconds
Alarm decibels:	65dB
Accuracy:	±4s in 24hrs
Battery:	1 x LR44 (included)
Operating conditions:	-10 to 50°C
CE & RoHS/WEEE:	Compliant
Hazard Information (SDS):	See www.brannan.co.uk for information
LCD dimensions:	12mm x 34mm
Product dimensions:	Net: 88mm x 48mm x 18mm Gross: 95mm x 53mm x 20mm
Product weight:	Net: 43.5g Box & instructions: 8.5g Gross: 45g



Description	Barcode	Product no
Slimline timer	5 011405 282132	28/213/0



View me on the website



brannan.co.uk

Measurement for life.

