

Brannan Product Data Sheet

PDS0292.20.REV3

Compact, Non-Contact, Infrared Body Thermometer



Non-contact and hygienic infrared thermometer designed to take the body temperature of a person regardless of room temperature. Designed to eliminate the chance of cross contamination, it is ideally suited for use in schools, businesses or other environments where a large number of people are coming into contact with one another.

Product Features:

- Body and object modes
- °C and °F temperature reading
- Last 10 readings memory storage
- Supplied in protective pouch, packed in individual box
- High temperature alarm
- 3 colour backlit LCD for easy reading
- Auto power off function

Product Specifications:

Effective Distance:	5cm to 8cm
Body Temperature Range:	34 to 43°C (93 to 109.4°F) <35°C / 95°F : "LO" displays <42.9°C / 109.2°F : "HI" displays
Body Temperature Accuracy:	+/-0.2°C / 0.4°F
Surface Temperature Range:	0 to 93°C (32 to 199.4°F) <0°C / 32°F : "LO" displays <93°C / 199.4°F : "HI" displays
Surface Temperature Accuracy:	+/-1°C / 1.8°F
Resolution:	0.1°C/°F
Memory:	Storage of last 10 readings
Measurement Time:	<1 second
Display:	3 colours backlit LCD (red, green, orange)
Operating Conditions:	16 to 40°C <85%RH (non-condensing) "ERR" displays when outside working conditions
Storage Conditions:	-20 to 55°C <85%RH (non-condensing)
Battery:	2 x AAA (included)
Power Consumption:	During power off: <10uW During measurement: <30mW
Power Level Indication:	Indication for low power level
Auto power off:	After 30 seconds of inactivity
CE & RoHS/WEEE:	Compliant
Hazard Information (SDS):	See www.brannan.co.uk for information
Product Dimensions:	Packaged: 165mm x 88mm x 48mm Unpackaged: 136mm x 76.5mm x 38mm
LCD Dimensions:	25mm x 23mm
Component Weights:	Net: 89g (including batteries) Protective Bag: 9g Display Box & Instructions: 33g
Gross Weight:	131g



Description	Barcode	Product No
Compact Infrared Body Thermometer	5 011405 114051	11/405/3

www.brannan.co.uk



Measurement for Life.

