

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

TRIPLE TIMER AND CLOCK

Three x 20 hour countdown timers with clock and memory function. Count up timer in 1/100 seconds. 12/24 hour clock. With stand and magnetic fixings.

Product features:

- 3 x countdown timers
- Count up timer and clock
- Large, clear LCD
- Desk stand and magnetic fixing.
- Complete with battery and instructions

Product specifications:

Timer count:	23 hours, 59 minutes, 59 seconds
Operating conditions:	-20 to 60°C < 90%RH (non-condensing)
Storage conditions:	-30 to 70°C < 75%RH
Battery:	1 x AAA (included)
CE & RoHS/WEEE:	Compliant
Hazard information (SDS):	See brannan.co.uk for information
LCD dimensions:	45 x 40mm
Product dimensions:	Gross: 185mm x 125mm x 38mm Net: 86mm x 56mm x 22mm
Product weights:	Net: 79g (including battery) Display card: 5g Display blister: 15.5g Gross: 100.5g



Description	Barcode	Product no
Triple timer and clock	5 011405 282071	28/207/0

