## **Product data sheet**





<u>+44 (0)1946 816600</u>



brannan.co.uk

## **FOUR WAY TIMER & CLOCK**

Four timers in one compact unit. 99 hour countdown and up. 12/24 hour clock. Counter automatically continues to count up when countdown time reached.

## **Product features:**

- Four timers in one compact unit
- Count up and down
- High, low and mute alarm volume
- Large, easy to read display
- Complete with hanging hook, magnetic back and stand



## **Product specifications:**

Timer keeping: Hour, minutes, seconds, AM/PM timer indicators (T1 to T4)

Timer count: 99 hours, 99 minutes

12/24 hour clock

+/- 30 seconds/month Time accuracy standard:

0 to 40°C < 95%RH (non-condensing) Operating conditions:

2 x AAA (included) **Battery:** 

CE & RoHS/WEEE: Compliant

**Hazard information (SDS):** See brannan.co.uk for information

LCD dimensions: 23mm x 55mm

**Product dimensions:** Net: 70mm x 88mm x 20mm

Gross: 100mm x 117mm x 23mm

**Product weight:** Net: 86g (including battery)

> Display box: 17g Gross: 102.5g



Description	Barcode	Product No
Four way timer and clock	5 011405 282200	28/220/0



