

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

ADVANCED DIGITAL VETERINARY THERMOMETER

Easy to use waterproof, electronic thermometer with a flexible probe and high speed response for quick and accurate reading.

Supplied with chart giving Hobday's specification (the normal body temperature for a variety of animals).

Product features:

- Compact, accurate and durable
- Choice of Celsius only or Celsius and Fahrenheit switchable readings
- Fast response, waterproof thermometer
- Non hazardous
- Complete with battery, instructions and Hobday chart



Product specifications

Temperature range:	32°C to 42.9°C or °C&F
Ingress protection:	IP22
Accuracy:	+/-0.1°C between 35.5°C and 42.9°C, otherwise +/-0.2°C
Resolution:	0.1°C
Approx response time:	60 seconds
Expected lifespan:	Approx 5 years
Battery:	1 x LR41 or equivalent (included)
Approx battery life:	100 hours continuous use (2+ years when idle)
Storage temperature:	-25°C to +55°C / <95% RH
Operating temperature:	+5°C to +35°C / <80% RH
CE & RoHS/WEEE:	Compliant
Hazard information (SDS):	See brannan.co.uk for information
LCD dimensions:	7mm x 20mm
Product dimensions:	Net: 133mm x 22mm x 12mm Gross: 167mm x 56mm x 18mm
Product weights:	Net: 19g (including batteries and probe case) Gross: 25g (including chart)

Hobday's Specification

Man	Cat	Horse	Dog	Cow	Pig	Sheep	Fowl
37°C	38°C	38°C	38.8°C	38.9°C	39.5°C	40°C	41.7°C

Description	Barcode	Product No
Advanced thermometer - °C only	5 011405 110350	11/035/2
Advanced thermometer - °C & F	5 011405 110350	11/035/3

View me on the website

