

# Product data sheet

✉ [sales@brannan.co.uk](mailto:sales@brannan.co.uk)

☎ +44 (0)1946 816600

🌐 [brannan.co.uk](http://brannan.co.uk)

## FOREHEAD STRIP TEMPERATURE INDICATOR

Quick and safe, forehead strip thermometer using heat sensitive liquid encased in plastic film. Ideal for small children or for adding into first aid kit.

### Product features:


- Simple to use and easy to read
- Quick and safe
- Supplied with storage sleeve
- Made in EU

### Product specifications:

Temperature range:	+37 to +40°C&F
Indicators:	Every 1.0°C/°F
Accuracy:	+/- 1.0°C
Storage & operating conditions:	16 to 40°C (60 to 104°F) at 15 to 95% RH*
CE & RoHS/WEEE:	Compliant
Product dimensions:	Strip: 89mm x 19mm Sleeve: 111mm x 79mm x 3mm
Product weight:	Strip: 0.5g Gross: 9g

### \*Note:

After use, store thermometer in the protective sleeve and keep away from Ultraviolet/sunlight, dampness, or heat. The thermometer can be cleaned by gently wiping with a soft cloth or tissue dipped in water. Do not immerse thermometer in liquids. Optimal storage and usage between 16 to 27°C (60 to 80°F) at 40 to 60% RH.



Description	Barcode	Product No
Forehead temperature indicator	5 011405 114662	11/466/9

