

THERMcert™ range

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

PRE-CALIBRATED MAX MIN, INDOOR & OUTDOOR THERMOMETER

High accuracy thermometer pre-calibrated at two fixed temperature points as part of the Brannan THERMcert™ range.

White moulded case with two displays for outdoor and indoor temperature. Each display has maximum and minimum recording facility and can be either Celsius or Fahrenheit. Complete with stand and external probe. Replacement probe also available.

Recalibration of your thermometer is recommended on an annual basis to ensure continued accuracy. Please contact our sales team for information.



Product functions:

- High accuracy thermometer supplied with traceable calibration certificate
- Measures indoor and outdoor temperatures
- Freezing alert
- 2m probe cable for remote reading
- Desk stand and hanging hook
- Battery and full instructions included

Product specifications:

Temperature range:	Indoor (internal reading): 0 to +50°C&F Outdoor (external probe reading): -50 to +70°C&F
Ambient operating temp.:	-10 to +50°C&F
Divisions:	0.1°C/F
Accuracy:	+/-2°C between -20 & 0°C +/-1°C between 0 & +50°C +/-2°C between +50 & +70°C
Probe:	Thermistor Probe IP67
Sampling rate:	10 seconds
Battery:	1 x AAA (supplied)
CE & RoHS/WEEE:	Compliant
Hazard information (SDS):	See brannan.co.uk for information
Cable length:	2500mm
LCD dimensions:	35mm x 52mm
Product dimensions:	Net: 98mm x 74mm x 17mm Gross: 275mm x 240mm x 34mm
Product weights:	Net: 79g Display box: 33.5g Gross: 128.5g (incl. certificate and packaging)



Description	Calibration points	Product no
THERMcert™ pre-calibrated max min indoor & outdoor thermometer	2 fixed points at 0 and 35°C	TC-019
Recalibration product reference: TC-019-RECALIBRATION		

