

# THERMcert™ range

## Product data sheet

✉ [sales@brannan.co.uk](mailto:sales@brannan.co.uk)

☎ +44 (0)1946 816600

🌐 [brannan.co.uk](http://brannan.co.uk)

### PRE-CALIBRATED SCOOP BATH THERMOMETER

Scoop bath thermometer pre-calibrated at two fixed temperature points as part of the Brannan **THERMcert™** range.  
Easy to use instrument to check bath water temperature by sampling.  
Dual scale with zones indicating recommended temperature.

#### Product features:

- Supplied with traceable calibration certificate
- Clear Celsius and Fahrenheit scale
- Zones for recommended bath temperatures
- Made in UK

#### Product Specifications:

<b>Range:</b>	-25 to +50°C&F*
<b>Accuracy:</b>	+/- 1.0°C @ ambient
<b>Divisions:</b>	1°C/2°F
<b>CE/ROHS/WEEE:</b>	N/A
<b>Hazard information (SDS):</b>	See <a href="http://brannan.co.uk">brannan.co.uk</a> for information
<b>Product material:</b>	White scale: high impact styrene Blue case: polypropylene Glass tube: kerosene filled
<b>Product dimensions:</b>	Net: 225mm x 40mm x 35mm Gross: 340mm x 250mm x 45mm
<b>Component weights:</b>	Plastic: 41.5g Glass: 6g
<b>Product weight:</b>	Net: 42g Display box weight: 15.5g Gross: 78g (incl. certificate & packaging)

\* Temperature ranges may have different upper and lower temperature extremes than stated due to manufacturing processes.



Description	Calibration points	Product no
<b>THERMcert™</b> pre-calibrated scoop bath thermometer	2 fixed points at 0 and 35°C	TC-014

**TRACEABLE CALIBRATION CERTIFICATE**

ISSUED BY: S. Brannan and Sons Limited

CALIBRATION DATE	CERTIFICATE NUMBER	SERIAL NUMBER

S. Brannan and Sons Ltd.  
Leconfield Industrial Estate  
Chichester Road  
Cumbria  
ENGLAND  
CA20 3SE

**Brannan**  
Thermometers  
and Gasps

Telephone: 01946 816600  
Fax: 01946 816609  
Email: [sales@brannan.co.uk](mailto:sales@brannan.co.uk)  
Web: [www.brannan.co.uk](http://www.brannan.co.uk)  
Registered No: 236123 ENGLAND

Instrument Type: SCOOP BATH THERMOMETER  
Range: -25/+50°C&F

Thermometer Reading	Calibration Temperature	Emergent Column Temperature	Standard Instrument References
0.0	0.0	N/A	246801/1 246801/2
34.9	35.0	N/A	246801/1 246801/2

Signed:  Calibration Assistant and Approved Signatory  
S. Brannan  
Laboratory Temperature 20°C +/- 5°C  
Laboratory Humidity 20% to 90%

S. Brannan and Sons operate an ISO9001 certified Quality Management System in order to ensure the accuracy and reliability of every instrument we produce.  
The readings were determined by comparison with traceable standard instruments. The reference numbers of which are shown above.  
We recommend that periodic recalibration should be carried out on this instrument yearly, but the recalibration frequency to be determined by the user depending upon the frequency and conditions of use, as well as the level of accuracy required by the process.