

THERMcert™ range

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

PRE-CALIBRATED WEATHERPROOF DIGITAL MAX MIN THERMOMETER

High accuracy thermometer pre-calibrated at two fixed temperature points as part of the Brannan **THERMcert™** range.

Robust, water resistant, digital max min thermometer with triple display simultaneously showing maximum, minimum and current temperature. Suitable for indoor or outdoor use. Hygienic and easy to clean – this is perfect for laboratories, catering and more.

Product functions:


- High accuracy thermometer supplied with traceable calibration certificate
- Brannan unique weatherproof design
- Maximum, current and minimum readings simultaneously
- High accuracy
- Clear, easy to read LCD
- Hanging hook

Product specifications:

- Temperature range:** -20°C to 50°C/°F (switchable)
Accuracy: +/-1°C
Ingress protection: IP66*
Display resolution: 0.1°C/0.1°F - current temperature
 1.0°C/2.0°F - max/min temperature
- Sampling rate:** Every 10 seconds
Calibration: 2 points @ 0 and 35°C
Battery: 1pc x AAA
CE & RoHS/WEEE: Compliant
LCD dimensions: 15.3mm x 52.3mm
 (current temperature display)
- Product dimensions:** Net: 210mm x 70mm x 20.5mm
Product weight: Net: 150g (including batteries)

Note: Current temperature LCD reads “LL.L” at -20°C and “HH.H” at 49.9°C

* IP rating determined as per BS EN 60529. Test duration: 3 minutes.



Description	Calibration points	Product no
THERMcert™ pre-calibrated weatherproof max min thermometer	2 fixed points at 0 and 35°C	TC-024
Recalibration product reference: TC-024-RECALIBRATION		



Measurement for life.

